

Dreamwire™

High Performance Guidewire



RX Biliary System™
A difference
you can feel.™

**Boston
Scientific**
Delivering what's next.™


Dreamwire™

High Performance Guidewire

DREAM ACCESS

Facilitates deep cannulation of targeted ducts


- 10cm hydrophilic Dream Tip™ End is designed for soft, smooth, access and excellent tactile feel
- Choice of angled or straight Dream Tip End
- Endoglide™ Coating promotes smooth tracking and enhances tactile sensation
- Forms an Alpha Loop to facilitate effective access into the targeted duct



DREAM EXCHANGE

Designed to improve procedural efficiency and easy RX Biliary device exchange

- 260cm and 450cm guidewire lengths available to meet physician preference for short-wire or traditional over-the-wire technique
- Kink-resistant Nitinol standard and stiff shaft, engineered for wire control and manipulation capabilities
- Striped movement markers designed to indicate wire positioning during wire advancement and device exchange
- Side loading Torque Device designed for effective guidewire tip manipulation (Angled-tip only)
- Guidewire may be left in place during sphincterotomy



ORDERING INFORMATION

| Product Number | Description | Length (cm) | Dream Tip Style | Shaft | Diameter (in) |
|----------------|------------------|-------------|-----------------|----------|---------------|
| M00556101 | Dreamwire ST | 260 | Straight | Standard | 0.035 |
| M00556111 | Dreamwire ANG | 260 | Angled | Standard | 0.035 |
| M00556121 | Dreamwire ST SS | 260 | Straight | Stiff | 0.035 |
| M00556131 | Dreamwire ANG SS | 260 | Angled | Stiff | 0.035 |
| M00556141 | Dreamwire ST | 450 | Straight | Standard | 0.035 |
| M00556151 | Dreamwire ANG | 450 | Angled | Standard | 0.035 |
| M00556161 | Dreamwire ST SS | 450 | Straight | Stiff | 0.035 |
| M00556171 | Dreamwire ANG SS | 450 | Angled | Stiff | 0.035 |

Packaged box/2

Indications, Contraindications, Warnings and Instructions for Use can be found in the product labeling supplied with each device.
CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

Boston Scientific

Delivering what's next.™

Boston Scientific Corporation
One Boston Scientific Place
Natick, MA 01760-1537
www.bostonscientific.com

Ordering Information
1.800.225.3226

© 2009 Boston Scientific Corporation
or its affiliates. All rights reserved.

DVG2420 2.5M April 2009